



Complete Machine QA



iba-dosimetry.com

## Your natural selection

### Platform based Machine QA

- ✓ Protocol based machine QA (including TG-142)
- ✓ Generic tests and customizable protocols
- ✓ Full coverage of tests with a flexible scheduling tool to manage your tasks, resources & time
- Comprehensive analysis, archiving, and reporting tools
- ✓ Interface to myQA Cockpit for quick and easy access to all QA results and trends

All results can be exported to reports, saved in the



benchmarked in the



and are displayed in the





# Your complete Machine QA solution

#### myQA Machines

is your complete protocol based machine QA on one platform. It offers a full coverage of tests related to dosimetry, safety, medical imaging, MLC QA... and more!

VMAT QA Plug-In



#### myQA Machines offers plug-ins for all of your QA needs:

#### **Dosimetry Plug-In**

- Perform automated dosimetry tests with the StarTrack\* or MatriXX detectors
- Acquire all key beam parameters in just one shot (dose output, profile analysis, energy verification)
- Analysis of main axis and diagonals (field size, symmetry, flatness, center, penumbra, light field)

# ibc StarTrackt

#### **MLC QA Plug-In**

- Perform automated MLC stripe tests ('picket fence test')
- ✓ EPID images & films can be analyzed to determine leaf position accuracy & MLC transmission characteristics
- ✓ Identifies whether any leaf is not in tolerance and which leaf number failed
- Every pixel of every strip is fitted to a modified Lorentzian curve

#### VMAT QA Plug-In

- Perform automated dynamic
  MLC QA and VMAT QA
- Automatic test analysis to verify accurate dose delivery using different dose rates, gantry speeds, and MLC leaf speeds
- Ensures that the changing rates and speeds during delivery do not adversely affect the delivered dose



#### **EPID QA Plug-In**

- Perform automated imaging QA for planar imaging (kV and MV)
- ✓ Fully automatic with all calculations performed in 5 seconds or less
- Compatible with all common imaging phantoms

#### **CBCT QA Plug-In**

- Perform automated imaging QA for CT and CBCT, including contrast, contrast to noise ratio, uniformity, HU deviation, spatial resolution, imaging scaling, and more!
- Compatible with all common imaging phantoms

#### PLUS



## Measurement Application for Star*Track*\* and MatriXX

- ✓ Instant display of results and real-time analysis (e.g.For beam steering, start-up behavior)
- Allows the Linac technician to setup unscheduled tests
- ✓ Import and export of measurements (ASCII)
- Profile comparison

*myQA* has given me the ability to have full control of my data by connecting all QA applications onto one platform and into one central database.

This platform allows me to cross-check data from within the various softwares, and offers the unique opportunity to connect with peers from around the world and benchmark QA data using the myQA Cloud. With myQA, the quality assurance becomes schedulable - in every sense of the word. For example, right now I am programming QA activities after Linac interventions using subgroups from periodic QA tasks, according to technical service recommendations. Another highlight for me is the myQA Cockpit. This web-based dashboard will allow us to quickly retrieve machine QA status updates anywhere in the department. myQA is truly an all-in-one solution.



Europe, Middle East, Africa | +49 9128 6070 North America, Latin America | +1 901 386 2242 Asia Pacific | +86 10 8080 9288 dosimetry-info@iba-group.com





iba-dosimetry.com