





# Your Global Quality Assurance Platform

All-in-One. All Connected. All Secure.

www.iba-dosimetry.com





✓ All QA applications integrated

✓ Simple and clear reporting

# Missing the connectivity in your QA?





- **X** No link between data, software programs & hardware devices?
- **X** Limited data compatibility between all your QA applications?
- **K** Lack of data comparability?
- **X** Difficulty in exchanging best practices and know-how with peers?

# Need a connected solution?



- ✓ All applications connected anytime & anywhere
- ✓ Connection to know-how & best practice
- ✓ Network with peers for the exchange of "tips & tricks"

# Accessibility & comparability of data, including benchmarks

# Looking for peace of mind in your QA?





- **X** Difficulty in following evolving technologies and protocols, and complying with regulations?
- **X** Unsure about best practice and most efficient QA procedures for each modality?
- **Worried about the safety of your data?**

# **Need a more supportive solution?**



*wy*QA<sup>™</sup> secure in all your QA aspects!



Secure Operations: Regular and automatic software updates to always remain up-to-date & compliant \*



Secure Outcome: Training support with the IBA Academy and clinical partners; 24/7 support from the Service Portal and Service Team



Secure Data: Integrity with consistent data storage and access from your central database or the cloud







## **Personal benefits.**

- ✓ Time saved with an easy and efficient workflow, a common look and feel, and simple functionalities
- ✓ **Full overview** with central database and "one click report"
- **Better results** based on simple comparison of all your data
- ✓ Gain knowledge and best practices from IBA peers
- ✓ **Peace of mind** based on Service Portal and Maintenance Contract program

# **Financial benefits.**

- ✓ Scalable solution: Simple tailoring to your specific needs
- **Regular automatic updates** protecting your investments with guaranteed costs for the contract period\*
- ✓ More patients treated with less QA effort
- ✓ Speed up your workflow with increased efficiency and connectivity
- ✓ Quick learning and implementation of new treatment modalities

# Strategic benefits.

- **Safer** patient treatments
- ✓ Integrate the whole department and satellite hospitals, and get connected with your peers!
- ✓ Global settings and common algorithms giving consistent results
- ✓ Maintenance contracts and training programs always keeping you up-to-date
- ✓ **Personalize** your QA platform exactly to your needs

**Treatment Safety.** 

Gain efficiency.

Be up-to-date.

Save time.

**Compare data**.

Make your QA personalized!

Treat more patients.

Exchange tips & tricks.

Save money.

Feel confident.













Get started now. Grow your Global QA platform with your needs!



# $\mathcal{W}_{\mathcal{V}}$ One Platform. Multiple Applications.



AT A GLANCE. Your brand new Global QA Platform integrates patient QA as well as machine QA and a Cockpit, including central user & equipment management. This platform enables you to access the different software modules and all data via one uniquely designed and intuitive interface, from different devices.

90% 🗸 05% A 05% × # Test Results ITrend Analysis Linac 1 🗶 Task Daily CA Cockpit Jask Weekly J Task Monthly panned J Task Yearly fidled / Linac 2 -IN COLUMN 2 IN COLUMN 2 passed / Linac 3 passed my QA Machines / Linac 4 -MatriXX Family passed -------Contract State State passed -----passed ( second . . my QA Cockpit Address, Contact of Female States, Sta Don't miss any key information 1))) 1111 Status **overview** of your clinic, your department, all your patients treatments plans and all your devices Star Track ✓ Traffic light style **Quick access** to your QA parameters and detailed test results my QA Accept **Track your data:** Trend Analysis and Statistics Browser application requiring **no local installation**, running on tablets and accessible from anywhere in your clinic 16 8 a 4 a **Blue Phantom** Family Just want to measure? Use the FAST TRACK !

Connect to your Star Track or MatriXX

✓ Start Measurement

✓ Instant results and real-time analysis

#### **Patient Plan Verification**

**Workflow optimized** plan verification in just 4 steps:

#### Import

**Patients** 

3

myQA









**3** Verify



✓ Increase your efficiency and reduce the patient QA time ✓ IMRT and rotational cases, global and local gamma Review the status of your patient verification in the Cockpit!

### Protocol based Machine QA

- ✓ Generic tests and customizable protocols (expected results, references & tolerances, frequency)
- ✓ Complete machine QA according to TG-142
- Dosimetry with Startrack or MatriXX
- CBCT QA, EPID QA
- MLC QA based on log files or images
- ✓ Scheduler: Manage your tasks and your time
- ✓ Import water phantom data to use as reference

✓ Machine QA in 3 steps:



## **Commissioning & annual Linac checks**

✓ Use the power of the Blue Phantom system for establishing vour high accuracy profiles

 Easily integrate your profiles e.g. into your *my* Amachines as a daily reference to your gold standard!



✓ Fastest, most Accurate, most Reliable





2







# **IBA** in a Nutshell

IBA offers innovative high-quality solutions and services with a focus on patient safety in cancer diagnosis and therapy.

#### Radiation Therapy: Fighting Cancer Safely and Precisely

- Industry-leading dosimetry and QA solutions that maximize efficiency and minimize errors for better outcomes
- ➤ Flexible soft tissue imaging markers (VISICOIL<sup>™</sup>) enable precise tumor targeting and tracking for pinpoint RT, IGRT and SBRT treatment accuracy and patient safety
- Undisputed leader in Proton Therapy delivery systems:
  Providing highly precise and effective radiation therapy

#### Medical Imaging: Safer Imaging, Earlier Cancer Detection

- Innovative Quality Assurance (QA) devices for x-ray dose and image quality checks, as well as QA of diagnostic displays
- > Patient dose monitoring solutions for x-ray imaging systems



#### dosimetry-info@iba-group.com

#### Europe, Middle East, Africa

IBA Dosimetry GmbH Bahnhofstr. 5 90592 Schwarzenbruck, Germany Tel.: +49 9128 607 0 Fax: +49 9128 607 10

#### North America, Latin America

IBA Dosimetry America 3150 Stage Post Drive, Suite 110 Bartlett, TN 38133, USA Tel.: +1 901 386 2242 Fax: +1 901 382 9453

#### Asia Pacific

IBA Dosimetry Asia Pacific No.6, Xing Guang Er Jie Beijing OPTO-mechatronics Industrial Park (OIP), Tongzhou District Beijing 101111, China Tel.: +86 10 8080 9288 Fax: +86 10 8080 9299







www.iba-dosimetry.com

myQA is pending release. Currently not available for sale in the USA